

uMED 20
Defibrillator Monitor

Saving lives through simplicity



7" TFT color screen



360J biphasic



3-lead ECG



50 mm recorder

Optional:



5-lead ECG/RESP



SpO₂



Pacer



AED

Technical Specification

Physical Specifications

Dimension 275 mm (w) x 205 mm (d) x 190 mm (h)
Weight ≤ 5.3 kg (basic config. with a battery and external paddles)

Display

Type TFT Color LCD
Dimensions 7 inch
Resolution 800 x 480 pixels

Power

AC Power

Line voltage 100 to 240 V~ (±10%)
Current 1.8 to 0.8 A
Frequency 50/60 Hz (±3 Hz)

Battery

Type 5600 mAh, rechargeable lithium ion battery pack
Capacity Monitoring mode: 6 hours, monitoring with 3-lead ECG
Defib mode: 300 times, 200 J discharge with three discharges per minute

Defibrillator

Waveform Biphasic truncated exponential waveform, with impedance compensation
Power on time Less than 2 seconds with fast startup mode
Charge time Less than 3 seconds to 200 J with a new, fully charged battery
Less than 7 seconds to 360 J with a new, fully charged battery
ECG recovery time Less than 2.5 seconds
Patient impedance range 25 to 300 Ω (external defibrillation)
Output energy 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 15, 20, 25, 30, 50, 70, 100, 120, 150, 170, 200, 300, 360J
Synchronized cardioversion Energy transfer begins within 60 ms of the QRS peak

ECG

Lead type 3-lead ECG
Lead selection I, II, III
Heart rate display Adult: 15 to 300 bpm
Pediatric: 15 to 350 bpm
Arrhythmia Yes
Alarms Yes
ECG size 2.5 mm/mV (x0.25), 5 mm/mV (x0.5), 10 mm/mV (x1), 20 mm/mV (x2), 40 mm/mV (x4), Auto

Sweep speed 6.25 mm/s, 12.5 mm/s, 25 mm/s, 50 mm/s
Patient isolation Type CF (defibrillation proof)

Data Storage

Patient profiles Max. 100 patients
Events Up to 500 events for one patient
Waveform storage Up to 10 hours of consecutive ECG waveform
Tabular trends 8 hours, resolution: 1 min
Data export Data can be exported to PC through USB flash memory

Environmental and Physical Requirements

Water/solids resistance IP44 (without external power)
Temperature Operating: 0 to 45 °C
Storage: -30 to 70 °C
Humidity Operating/storage: 15 to 95 % (non-condensing)
Altitude Operating/storage: -381 m to +4575m
Drop protection 0.75 m

www.mindray.com

P/N:ENG-uMED 20-210285X6P-20230118
©2023 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd. All rights reserved.

mindray
healthcare within reach



Simple, but Professional

As a monitor uMED helps recognize and prevent cardiac arrest from occurring; as a defibrillator uMED supports users through their individual level of proficiency in both manual and AED mode.

Manual Defibrillation with Clear 1-2-3 Steps




1. Select Energy

2. Charge


3. Shock

Intelligent ResQNavTM


With innovative ResQNavTM technology, uMED 20 can evaluate the proficiency level of clinicians, perform automatic rhythm analysis and intelligently navigate the rescue process step by step, even in manual defibrillation mode.



Paddle placement



Auto-analysis



Charging guide

External Paddles with Patient Contact Indicator

Buttons for energy selection, charging and shock delivery improve usability for clinicians.

Patient contact indicators both on the paddle and screen give visual status of paddle contact to ensure discharge efficiency.



Patient contact is good



Patient contact is not good




Patient contact is very poor



Quick Guidance, Easy Maintenance

Quick Guidance, Interactive Experience

uMED provides quick operation guide. With interactive immersive experience, the quick demo guide helps medical staff quickly grasp the key operations of the device.



Demo guide for key operations

Easy Maintenance, all at a glance

A defibrillator is on standby status for more than 95% of the complete life cycle. However, a lot of defibrillators on the market do not have comprehensive auto-test function and rely on manual check. How can device failures be detected and resolved in time? uMED 20 makes routine maintenance simple and safe, helping to save manpower, improve efficiency and ensure that the defibrillator is ready for use at any time.

- Device status is clear at a glance with comprehensive auto-test, no need to daily manual check.
- Quick access makes it easy for medical staff to view the whole test summary.
- Customized auto-test report helps the users only focus on what they really care about.
- Device failures can be simply resolved by users themselves, thanks to graphical visual troubleshooting guide.



Quick access to test summary



Reliable & Effective

Reliable & Durable Quality

To be reliable in facing a variety of possible accident, uMED 20 has passed strict safety and reliability tests. It is extremely durable and has a long life span.

- IP44 water-/dust-proof
- 6-surface 0.75m drop test without any additional protection
- Working temperature is 0~45°C, unaffected by extremes
- Large capacity battery, longer battery life, more durable



Fast & Effective Shock

uMED 20 is also equipped with 360J biphasic and QShockTM technologies, which further improve the effectiveness of defibrillation

- 360J biphasic technology with auto-compensation according to patient impedance, working for a wider range of people
- QShockTM technology brings an extremely fast defibrillation experience, for every second counts.

